

Averis Mini Centrifuge

Features include:

- Ideal for microfiltration and quick spin-downs
- Easy-to-use with dual start/stop function
- 5000 rpm fixed speed. Max RCF, 1360xg
- Versatile rotor compatibility with 0.2 to 2.0mL tubes and 0.2mL PCR strips up to 32 wells
- Whisper Quiet: 45dB
- Easy Change Rotor
- Accelerates and brakes in seconds
- Available in blue



Technical Specifications:

Max Speed 5000rpm
 Max. RCF 1360xg
 Rotor Capacity 0.2/0.5/1.5/2.0mL x 8;
 0.2mL x 32 PCR strips;
 or 0.2mL x 4 PCR 8 strips
 Run Time Continuous operation
 Fuse 0.5A, 250V (F)
 Power 110V, 60Hz, 130mA or 240V,
 50/60Hz, 65mA
 Noise Level 45dB
 Dimensions (D x W x H) . . . 160 x 170 x 122mm
 Weight 0.5kg
 Certification CE, cTOVus, FCC, MCA

Ordering Information:

CE ISO 9001

PART#	DESCRIPTION	QTY
AGV-D1008E	Averis Mini Centrifuge	1 Each

Exclusively Distributed by Astraglass Innovations

www.astraglassinnov.com • +1 (800) 955-0376 • sales@astraglassinnov.com