

## Averis Mini Dry Bath Incubator

Averis Mini Dry Bath Incubator is a reliable choice for precise temperature control. This compact dry bath offers efficiency and convenience, making it ideal for a variety of applications including enzyme reactions, DNA denaturation, cell incubation, and nucleic acid hybridization.

### Key Features & Benefits:

- **LCD display with simultaneous readout of time and temperature**
- **Small Footprint, 150 x 110 x 100mm — fits any crowded benchtop**
- **Simple touchpad digital control**
- **Accommodates 7 different aluminium block configurations**
- **Digital temperature control for precise, repeatable adjustments**
- **Temperature Deflection Adjustment**
- **Audible alarm on cycle completion**
- **12 V DC power supply**



*Blocks Not Included*

### Technical Specifications:

Power Supply . . . . .	DC 12 V
Power . . . . .	35 W
Temperature Range. . . . .	RT+5°C – 80°C
Time Range. . . . .	1 sec – 999 sec or 1 min – 999 min
Temperature Accuracy . . . . .	±0.5°C
Display Accuracy . . . . .	±0.1°C
Heating Time (25°C to Max.) . . .	≤12min
Ambient Temperature . . . . .	5°C ~ 35°C
Dimensions (WxDxH) . . . . .	150 x 110 x 100 mm
Weight . . . . .	0.5 kg

### Ordering Information:



PART#	DESCRIPTION	QTY
AGV-AS-01130-00	Averis Mini Dry Bath Incubator (Blocks not included)	1 Each
AGV-AS-01131-01	40 x 0.2mL Block	1 Each
AGV-AS-01131-02	24 x 0.5mL Block	1 Each
AGV-AS-01131-03	15 x 1.5mL Block	1 Each
AGV-AS-01131-04	15 x 2.0mL Block	1 Each
AGV-AS-01131-05	8 x 12.5mL Square Cuvette Block	1 Each
AGV-AS-01131-06	4 x 15mL Block	1 Each
AGV-AS-01131-07	2 x 50mL Block	1 Each

Exclusively Distributed by Astraglass Innovations

www.astraglassinnov.com • +1 (800) 955-0376 • sales@astraglassinnov.com