

### Averis Mini Orbital Shaker

**Features include:**

- Variable speed 50–250 RPM
- Small foot print
- Quiet operation
- Four platform options available
- Digital LED Display
- Timer 0–99:59 (mm:ss)



**Technical Specifications:**

Platform.....	23 x 27cm
Speed .....	50 to 250 rpm
Timer.....	0–99:59 (mm:ss)
Unit Dimensions (w x h x d) .....	28 x 11 x 27 cm
Universal Voltage Input.....	100–240V, 50/60HZ
Safety Certifications .....	EN61010-1:2001, CE

**Each Model Includes:**

- Shaker
- Table top
- Screw Driver and 4 screws
- Power cord (115V or 230V depending on catalog number)
- User Manual

**Ordering Information:**



PART#	DESCRIPTION	QTY
AGV-AS-04030-00	Averis Mini Orbital Shaker	1 Each
ACCESSORIES & REPLACEMENT PARTS		
AGV-AS-04021-01	Platform for 12 x 100 to 150mL Flasks	1 Each
AGV-AS-04021-02	Platform for 6 x 250 to 300mL Flasks	1 Each
AGV-AS-04021-03	Universal Platform with Adjustable Bars	1 Each
AGV-AS-04021-04	23 x 27cm Flat Table	1 Each