

Averis Mini Stir Plate

Averis Mini Stir Plate features a sleek, lightweight and compact design. The built-in micro motor provides smooth, reliable stirring. A 1-year limited warranty is included with every unit. Single-position unit with a maximum capacity of 1 litre.

Features include:

- **Small Footprint:** 116 x 116 x 35mm preserves valuable benchtop space
- **Recommended Stir Bar Size:** 8mm x 30mm
- **Diameter of Plate:** 116mm
- **Maximum Mixing Capacity** 1 litre
- **Adjustable Stirring Speed:** 200–2000 rpm, continuously variable via dial control
- **Control Mode:** Single Knob
- **Direction of Rotation:** Clockwise
- **Magnon Size Range:** 20–30mm
- **Motor:** A 5W Micro Motor provides consistent stirring, offers smooth, quiet and efficient operation
- **Weight:** 460gm
- **Warranty:**
Averis offers a 1-year limited warranty



Ordering Information:

CE ISO 9001

PART#	DESCRIPTION	QTY
AGV-FSF-MN	Averis Mini Stir Plate	1 Each