

Averis Pipette Controller

Features include:

- Averis Pipette Controller enables single-handed operation with minimum effort
- Compatible with most plastic or glass pipettes from volumes of 0.1–100mL
- Easily adjustable speed control
- 0.45µm replaceable hydrophobic filter
- Lightweight with reliable performance
- High capacity Li-ion battery for extended operation time

Technical Specifications:

Aspiration Speed	25 mL in < 7 seconds
Dispensing Speed.	Motor / Gravity
Battery.	Lithium-ion
Battery Service Life.	8 Hours of Intermittent Use
Charging Time.	2–3 Hours
Pipette Types	Glass or Plastic Pipette (0.1–100mL), Pasteur Pipettes
Filter.	0.45µm Hydrophobic (replaceable)
Weight	200g



Ordering Information:

CE ISO 9001

PART#	DESCRIPTION	QTY
AGV-Levo Me	Averis Pipette Controller, 110V	1 Each

Exclusively Distributed by Astraglass Innovations

www.astraglassinnov.com • +1 (800) 955-0376 • sales@astraglassinnov.com